Feasibility of PT, aPTT, and Fibrinogen on a Multifunctional Digital Microfluidic Cartridge

BACKGROUND

- Measurement of multiple coagulation parameters quickly is vital to correcting coagulopathy and improving outcomes in severe trauma patients.
- There are no rapid, near-patient technologies that can perform multiple coagulation assays simultaneously on a single consumable.
- We are developing a rapid (<15 minutes), point-of-care multifunctional digital microfluidic (DMF) platform to test critical coagulation markers prothrombin time (PT), activated partial thromboplastin time (aPTT), and fibrinogen – using low volume samples (<50 μ L).

METHODS

- 10 µL samples are loaded onto a DMF cartridge and then dispensed to combine with PT, aPTT, or fibrinogen reagents.
- Reaction droplets are then shuttled within the cartridge where the droplet motion is monitored through electrical impedance to measure the rate of change in viscosity of each reaction. Droplet viscosity changes at a rate proportional to the amount of coagulation precursors present, allowing for the calculation of clotting time and activity.
- To determine accuracy of the PT and aPTT assays, five proficiency testing plasma samples (College of American Pathologists [CAP], Northfield, IL) with known values for normal, mid-level abnormal, and high-level abnormal clotting were tested on the DMF platform.
- For method comparison studies, contrived plasma samples containing K2 EDTA and unfractionated heparin to simulate prolonged time were tested on the DMF platform for PT and aPTT assays respectively and compared to results from an FDA-cleared comparator method. CAP and patient plasma samples (LabCorp, NC) were tested on the DMF platform and an FDAcleared comparator for the fibrinogen assay.



Figure 1. Baebies' point-of-care platform for rapid coagulation testing

Pamula VK, Ceasrine A, Ruddiman C, Caudle SB, Ng R Baebies, Inc., Durham, NC 27703



RESULTS





Figure 3. Method comparison studies against FDA-cleared comparators using A) contrived plasma for PT (n=37), B) contrived plasma for aPTT (n=7), and C) CAP and patient plasma samples for fibrinogen (n=26). The DMF PT, aPTT, and fibrinogen assays correlated to the comparator method with R² values of 0.80, 0.95, and 0.96, respectively.



Research reported here was supported by the National Heart, Lung and Blood Institute (NHLBI) of the National Institutes of Health under award #R44HL169045. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

These assays are not available at this time for sale or use in any territory

PT, aPTT, and fibrinogen results correlated well with comparator testing in plasma samples, and we are currently adapting the assays for use on whole blood

DMF coagulation testing with a rapid turnaround time can aid in guiding hemostatic resuscitation in severe trauma patients.

haphips